

STATdx

Quick, Confident Diagnosis.

Jasmine Tsai 蔡雅莉 Senior Account Manager j.tsai@elsevier.com



ELSEVIER


STATc

New Content in STATdx

Our content has been recently updated. Be on the lookout for:

- Updated content in Breast including:
 - Over half of the 4,000+ images have been updated, including imaging findings complemented by histopathologic and clinical correlates of the spectrum of breast cancer
 - New topics such as gender reassignment, expanded information on disease-causing mutations and updated risk assessment, and in-depth coverage of malignancy in pregnancy and male breast disease, providing a helpful resource for such less common but important topics
 - Multi-perspective analysis of updated breast cancer screening guidelines, including but not limited to United States Preventive Services Taskforce, American Cancer Society, and American College of Radiology guidelines
- New Brain Dx topics written for Zika Virus Infection, Metronidazole Toxicity, and Gadolinium Deposition
- CT and MR in Cardiology has been updated with 155 new or revised topics

GOT IT



Feature Highlights—Rely on the Experts

STATdx, written by renowned radiologists in each specialty, provides comprehensive decision support you can rely on, including:

- □ Over 200,000 image examples of X-ray, CT, MRI, and ultrasound
- □ Over **4,300** common and complex **diagnoses**
- ☐ More than **1,300 differential diagnosis** modules
- □ 300 comprehensive normal imaging anatomy
- 18,000+ easily-sortable patient cases
- ☐ 195 basic and advanced intervention procedures
- □ Statdx(影像資料庫) covid-19 專區 https://player.vimeo.com/video/405035237
- □ The Power of STATdx and ClinicalKey Anne Osborn, MD, FACR https://www.statdx.com/video/



Feature Highlights

1.RADTools (RADTools)

醫事放射師常用工具之彙整,包含TNM和癌症分期表等圖表,各式診斷程序, 對齊角度之參考與分類資訊,以及各種計算工具。

2.主題預覽 (Topic Previews)

在目錄、搜尋結果、以及診斷模組中瀏覽主題標題時,可預覽相關內容。

3. 隨時比較 (Compare Anywhere)

可並列對照兩個診斷影像以上,方便您快速而清楚的交互參考。

4.簡報製作


每張圖片均可自動匯出成PPT,並附上圖片說明,讓您輕鬆準備演講題材。

5.關鍵字 link out 搜尋關鍵字可串連至 Pubmed/radiology pdf/google


6.病例應用 含括18000筆的病例,可以深入探討應用

7.專家建議 可快覽專家建議,方便臨床應用

AMIRSYS STATdx 影像資料庫簡介



Search by keyword



Google Scholar

PDF Search

Compare

比對不同的疾病診斷圖像




liver scar





Filter by Category





ddx Hepatic Mass with Central Scar

Ultrasound by Gregory E. Antonio, MD, DRANZCR, FHKCR

12 images Updated 8/4/14



Filter b

Clid

to pr

... necrosis, and calcification, which may simulate central scar; Color Doppler shows hypervascular tumor supplied by hepatic artery. Hepatic



ddXLiver Mass with Central or Eccentric Scar

Gastrointestinal by Michael P. Federle, MD, FACR





... Nearly isodense (isointense) to liver on all ... Hepatic Cavernous Hemangioma: Large hemangiomas (> 5 cm) commonly have fibrotic scar that may ...



dX Focal Nodular Hyperplasia

Gastrointestinal by Michael P. Federle, MD, FACR

29 images: 19 references Updated 9/20/12





... or isodense to normal liver; Large draining veins -- hepatic veins. Delayed scans: Mass: -- isodense to normal liver; Central scar: Hyperdense (due to ...



dX Focal Nodular Hyperplasia

Ultrasound by Gregory E. Antonio, MD, DRANZCR, FHKCR

22 images : 3 references





... there is no significant mass effect or bulge in liver contour; ... Mass effect: Displacement of normal hepatic vessels and ducts; Central scar: Seen in ...



dX Fibrolamellar HCC

Gastrointestinal | by Michael P. Federle, MD, FACR

25 images: 11 references Updated 7/23/12





... Conventional HCC in noncirrhotic liver mimics FLC. Hepatic Cavernous Hemangioma. ... mass; Central decreased attenuation (scar), rarely with ...






5 images: 4 references





Compare

比對不同的疾病診斷圖像






RAD Tools-Calculators

	Calc	Calculators			
Î	œ.	Ankle Fractures Calculator	1 image : 1 reference Updated 04/16/15	୯	
	₩.	<i>▶</i> Bone Age Calculator	1 image : 2 references Updated 05/11/15	୯	
	ΞΔΞ	Bone Tumors Calculator	1 image Updated 05/11/15	୯	
		Intracranial Cysts Calculator	1 image Updated 05/11/15	C ₂	
1	*	Lung Cancer Staging Calculator	1 image : 1 reference Updated 05/11/15	C ₂	
	A	Pneumothorax Calculator	1 image : 2 references Updated 05/11/15	Co.	
1	550	Radiation Dose Calculator	1 image : 1 reference Updated 05/11/15	C ₂	
1		Renal Insufficiency Calculator	1 image : 1 reference Updated 05/14/15	C ₂	
1		Solitary Pulmonary Nodule: Bayesian Method	1 image : 2 references Updated 05/14/15	C ₂	
Î	Nachtables Seden Ag Notice Sarate Nationage Step	Solitary Pulmonary Nodule: Gould Method	1 image : 1 reference Updated 05/14/15	୯	
Î		Tumor Doubling Time: Diameter	1 image : 1 reference Updated 05/14/15	୯	
Î	3	Tumor Doubling Time: Volume	1 image : 1 reference Updated 05/14/15	୯	

RAD Tools-Tables


- AJCC Tables
- Neuro Tools
- Obstetrics Tools
- Chest Tools
- GI/GU Tools
- MSK Tools
- Pediatrics Tools

T1 Hyperintense Basal Ganglia



Key Facts

摘錄自AMIRSYS電子書系列的完整內容



Script



Parenchymal Brain Tumors

Paula J. Woodward, MD


Hide Images

Selected Images






Sagittal graphic shows a heterogeneous pineal region teratoma. There are cystic A and solid areas within the mass. Calcifications -> are the most specific sign of teratoma but are not always present. Most fetal brain tumors begin in the pineal region but grow so large that the point of origin is often not discernible.



Sagittal T2WI MR of a 3rd trimester fetus with a teratoma shows the mass compressing the cerebrum and stretching the brainstem and




Transverse ultrasound of a fetal brain shows a large, heterogeneous mass within the cranial vault completely destroying normal anatomic landmarks. Measurements showed marked macrocephaly.




Postmortem coronal T1WI shows complete replacement of brain tissue by a complex mixed signal intensity mass. Immature teratomawith primitive neural ectodermal tissue, cartilage, bone, intestinal mucosa, smooth muscle, and hemorrhage a was identified at autopsy.




Transverse ultrasound of a fetal brain shows a mixed cystic and



Gross pathology in the same case. shows a variegated, lobular mass with marked thinning of the rampining carehral ticcos




Axial T2WI MR of a 3rd trimester fetus shows a slightly hypointense,




Sagittal T2WI MR on day 1 of life shows a large heterogeneous mass with high signal cystic areas

圖像自動轉存為PPT



Bone Age Calculator







Elsevier's STATdx

- A diagnostic decision support system for radiologists



<u>USE CASE</u> Lyme's Disease

Outline

- Lyme's Disease
 - Clinical History
 - MRI Findings
 - Navigating STATdx
 - Primary Differential Diagnosis
 - Reference Images
 - About STATdx

Clinical History


- 35-year-old woman complaints of weakness on one-side of the face after getting back from a camping trip
- Complaints also include severe headache, neck stiffness, fever, chills, muscle/joint pain and fatigue
- Patient presents with multiple erythema migrans (skin rashes), round, outwardly expanding rash ("bull's-eye")

 Patient's relative also complaints that this woman is having concentration difficulties and unable to remember certain information lately


Confirmation of diagnosis requires: ELISA, PCR, MRI

MRI Findings.

- Axial FLAIR image demonstrates multiple foci of signal abnormality involving the periventricular white matter bilaterally.
- Contrast-enhanced axial T1WI does not demonstrate any enhancement.
- Contrast-enhanced, fat-suppressed axial T1WI demonstrates enhancement in the labyrinthine and anterior tympanic segments of the left facial nerve..








Search function helps user find relevant content

Users may search by documents or images





ddx Enhancing Cranial Nerve(s)


ddX Ring-Enhancing Lesions

Brain by Anne G. Osborn, MD, FACR

ddx Ring-Enhancing Lesions, Multiple

Pediatrics by Bernadette L. Koch, MD

... Aneurysm (Thrombosed); Other Infections: Tuberculosis; Fungal Diseases; Acquired Toxoplasmosis; Lyme Disease. Other Neoplasms:



Click the "Preview" i to preview an item's c

? Help/Support

CME


Primary Differential Diagnosis

Demyelinating disease

Vasculitis

Sarcoidosis

Chronic fatigue syndrome




Decision support s

 Axial FLAIR image demonstrates multiple foci of signal abnormality involving the periventricular white matter bilaterally,.

Contrast-enhanced axial Table
demonstrate
Unique imaging clues help radiologists choose which
Unique imaging clues help radiologists choose which
diagnoses to consider and compare

a suppressed a sup

Diagnosis: Lyme's Disease






Each diagnosis topic includes patient cases:

Each case includes demographics, history, case description, and the author/contributor.

The Description includes the author notes and annotations for the case.

All cases include numerous labeled images.



About STATdx

A point-of-care diagnostic decision support system for working and studying radiologists.

STATdx increases speed, accuracy and confidence in diagnosing complex imaging cases

More than 4,000 diagnoses written by the world's leading experts in radiology

STATdx includes over 1,000,000 images, including x-ray, CT, MR and ultrasound images

Nearly 20,000 supporting individual patient cases, authored by world-renowned experts in radiology

Unique 'Compare' feature allows a side-by-side comparison of up to three diagnoses simultaneously

